

8821Q-R USB FLASH DRIVE

The 8821Q-R USB Flash Drive contains the following items;

- Q-Lab Software - This software installer is located in the **Q-Lab Software** folder.
- Q-Lab Operation Manual (P/N 0010268001) - This Operation Manual is located in the **Operation Manuals** folder.
- Q-Lab Operation Manual (P/N 0010267000) - This Operation Manual is located in the **Operation Manuals** folder.

To access the files located on this 8821Q-R USB Flash Drive, insert the device into an available USB port on a PC or Laptop Computer and copy the folders and their contents to a convenient location on your local hard drive.

To install the Q-Lab Software, double-click on the file titled “Q-Lab.ms” that is located in the **Q-Lab Software** folder.

For more information on the installation and setup of the Q-Lab Software, see the included Q-Lab Software Operation Manual located in the **Operation Manuals** folder.

If you need assistance with this software, please
call: (800) 344-2412 or (317) 895-3600
P/N 0010268002 rev. 8/08

